

c-Myc-tag control-lysate

Description

The c-Myc-tag control-lysate was generated from transfected HEK293 cells, expressing a c-Myc tagged protein. The c-Myc-tag corresponds to a decapeptide (EQKLISEEDL) of the c-Myc proto-oncoprotein and belongs to the most widely used epitope-tags.

To generate the c-Myc-tag control-lysate HEK293 cells were transiently transfected with a mammalian expression vector coding for a cytoplasmic protein harbouring an amino-terminal c-Myc-tag. A whole cell lysate was prepared from these transiently transfected cells. Besides all other cellular proteins, this whole cell lysate contains the c-Myc-tag protein, which shows an apparent molecular weight of ~21 kDa upon SDS-PAGE.

In addition to the c-Myc-tag control-lysate, a mock-lysate was prepared from HEK293 cells, which were transfected with the empty expression vector.

Product

Whole cell lysate:

50 µl

Cell line:

HEK293 cells

Tagged protein:

Human cytoplasmic protein with an amino-terminal c-Myc-tag.

MW:

~21 kDa

Lysis buffer:

60 mM Tris/HCL pH 6.8, 2% SDS, 5% mercaptoethanol, 50 µg/ml bromophenolblue, 10% glycerol

Application:

For SDS-PAGE/Western Blotting 20 µl of the lysates should be applied per lane of a 15% polyacrylamide-gel.

Shipping and storage

Shipping:


The lysate is shipped in cold case.

Storage:

Lysate is stable for 1 week at 4°C.

For prolonged storage, the lysate should be stored at -80°C. Aliquot to avoid repeated freezing and thawing. At -80°C, the product is stable for at least 1 year from shipment.


Application in quality control



Control cell lysates
anti-c-Myc-tag blot

Western Blotting using the c-Myc-tag control-lysate and the mock-lysate.

HEK293 cells were transiently transfected either with the empty expression vector (mock-lysate; mock lys.) or with an expression vector encompassing the cDNA of a c-Myc-tag protein (c-Myc-tag control-lysate). 25 µl of the mock-lysate and the indicated amounts of the c-Myc-tag control-lysate (25, 20, 15, 10, or 0 µl) were separated by SDS-PAGE using a 15% polyacrylamide-gel and then transferred onto a PVDF membrane. The membrane was probed with the monoclonal anti-c-Myc-tag antibody from

 tag-tools (clone TT5; dilution 1:2,000).

The c-Myc-tag bound primary antibody was detected by incubation with HRP-coupled protein G (dilution 1:10,000) and visualized using chemiluminescence.

The X-ray film was exposed for 30 seconds. Whereas a specific ~21 kDa band is detected in the c-Myc-tag control-lysate at all volumes (10 - 25 µl), no binding of the anti-c-Myc-tag antibody to cellular proteins can be detected in the lane, where 25 µl of the mock-lysate was applied.

Use

The c-Myc-tag control-lysate and the mock-lysate are for quality control in research applications only. Not for diagnostic or therapeutic purpose.